

## RUMEX GLAUCOMA SHUNT, NON-VALVED

- Large surface area allows aqueous filtration and pressure relief by controlling long-term intraocular pressure (IOP), minimizing postoperative complications and medication needed in case of use of small plate design implants.
  - The design of the implant allows implantation of the large surface plate avoiding the risk of trauma, ensures secure plate position and fast healing.
  - Unique fenestration design allows the growth of fibrotic tissue through the fenestrations to help control bleb height and volume and minimize ocular motility disturbance.
- 
- Material: Medical grade silicone



**GlaucoDrain**  
Aqueous drainage device  
Антиглаукомное дренажное устройство

